

# A Cysteine-containing Tag and Complimentary Adapter Protein for Loading Contrast Agents Onto Targeting Proteins

Marina V. Backer, Vimal Patel, Brian T. Jehning and Joseph M. Backer

SibTech, Inc., Newington, CT 06111

## Abstract

Targeted imaging of disease-associated cell surface markers can lead to development of personalized treatment regiments. The lack of a reliable technology for "loading" contrast agents onto targeting proteins inhibits the advances in this area. Currently, loading is achieved by random chemical cross-linking of cargo to targeting proteins, which damages proteins, requires expensive customized development, and yields heterogeneous products.


We have developed a new platform technology for loading contrast agents onto targeting proteins that circumvents these problems. The strategy is based on expression of a targeting protein with a 15-aa humanized tag containing cysteine in position 4 (C-tag). The tag is fused to a targeting protein via a short G<sub>4</sub>S linker.

In our experience, C-tagged recombinant proteins produced in bacteria and refolded in typical red-ox buffers contain varying numbers of free thiol group (5-40%), while the majority of the C4 thiol groups requires "deprotection" under mild reducing conditions. After deprotection, C4 thiol group can be conjugated either directly to contrast agents (e.g. near-infrared dyes) or to contrast agent carriers (e.g. chelators for radionuclides, liposomes).


Alternatively, contrast agents or their carriers can be conjugated to a ~10 kDa humanized adapter protein that can bind to C4-tag. To avoid dissociation of the complex, adapter protein is engineered to form an intermolecular disulfide bond with C4 thiol group upon the complex formation. Importantly, formation of this disulfide bond does not require prior "deprotection" of C4 thiol group. We expect that these facile procedures for loading contrast agents onto proteins will provide new opportunities for targeted imaging.

## 1. Site-specific Cross-linking of Contrast Agents


Random crosslinking damages proteins






Fusion Tag for site-specific modification



Complimentary Adapter protein



## 2. Fusion Tag and Adapter Protein Are Parts of Human RNase I




Human RNase I can be reconstituted from two fragments: a 15-aa Tag and a 107-aa Adapter

Adapter protein is released by limited protease digestion of chimeric BH-RNase

High affinity interaction of Tag and Adapter is used for purification of Adapter from protease digestion mixture

## 3. Site-specific Modification of C-tagged Targeting Proteins



## 4. Payloads Used for Cross-linking


Chelators:  
HYNIC, DOTA  
Near-infrared dyes  
PAMAM dendrimers  
Pegylated phospholipids  
20 kDa and 40 kDa PEG

## 5. C-tagged Targeting Proteins

Vascular endothelial growth factor (VEGF)  
Proprietary single-chain VEGF  
Annexin V  
Catalytically inactive fragment of anthrax lethal factor (LFn)  
scFv single-chain antibody fragments

## 6. C-tagged VEGF for Imaging Tumor Neovasculature


Engineering C-tagged dimeric and single chain VEGF



**C-VEGF<sub>121</sub>** : human VEGF<sub>121</sub> was cloned into C-tag vector for bacterial expression.


**Single chain C-VEGF** (scVEGF): two 3-112aa fragments of human VEGF<sub>121</sub> were cloned head-to-tail into C-tag vector for bacterial expression

C-VEGF binds Adapter via S-S bond




After a 16-hrs incubation at 4 °C, samples of Adapter protein, C-VEGF<sub>121</sub> and their mixture at molar ratio of 5:1 were analyzed by SDS-PAGE under reducing or non-reducing conditions.

Various VEGF conjugates bind and activate VEGFR-2




Various VEGF conjugates rescue cells from VEGF-Toxin




**C-tagged VEGF driven *in vivo* fluorescent and SPECT imaging of 4T1 tumors is presented on posters # 303 and 585**

## 7. Conjugation of Liposomes to Targeting Proteins: Directly on C-tag or via Adapter. The story of C-tagged LFn.


Doxil liposomes can be bound directly to C-tag or via Adapter




LFn-Adapter conjugate forms within minutes of incubation




Only LFn-Ad-Doxil effectively kills RAW cells in the presence of PA (72-hr assay)



C-tagged LFn and LFn-Adapter are functionally active in protection of RAW cells from PA-mediated LFn toxicity



LFn-Doxil and LFn-Ad-Doxil are functionally active (2-hr assay)




RAW cells were plated on 96-well plates (10,000 cell/well). Twenty hrs later, targeted Doxil was serially diluted in complete culture medium with or without PA and added to cells in triplicate wells. After a 5-hr incubation at 37°C, Doxil-containing media were aspirated and cells were shifted to complete culture medium. Cells were quantitated after 72 hrs of total time by an MTT-based assay.

RAW cells were plated on 96-well plates (10,000 cell/well). Twenty hrs later, free or modified LFn was serially diluted in complete culture medium with indicated amounts of LF and PA, and added to cells in triplicate wells. After a 2-hr incubation at 37°C, cells were quantitated.


## 8. C-tagged Annexin V for Imaging Apoptosis

C-tagged annexin binds Adapter protein



Annexin-Adapter protein conjugate was detected in annexin/Adapter protein mixture after a 16-h incubation at 4°C. The conjugate was purified from unreacted components by ion-exchange chromatography (SP column) followed by affinity chromatography (Tag-bound column).

Annexin-Ad conjugate displays unaffected functional activity



Functional activity of annexin-Adapter was estimated by its ability to compete with FITC-annexin for binding to phosphatidylserine-displaying erythrocytes of stabilized human blood. Annexin binds to erythrocytes in the presence of Ca<sup>2+</sup> and is released in the presence of EDTA. Both annexin-Adapter conjugate and C-tagged annexin displaced FITC-annexin in a dose-dependent manner with IC<sub>50</sub> of 11±3 nM. Correct size recombinant human annexin V displays similar activity in this assay with IC<sub>50</sub> of 9±4 nM

## Conclusions

- A platform technology for standardized loading of imaging and therapeutic agents onto targeting protein is developed.
- The technology is based on a 15-aa cysteine-containing tag and a complimentary adapter protein and offers two strategies for conjugating payloads to targeting proteins.
- The technology works with different recombinant proteins that are expressed with the cysteine-containing tag.
- The technology works with various contrast agents and their carriers such as chelators, dendrimers, liposomes.